

CS4

FURNACE APPARATUS
FOR KLY5-A AND MFK1-FA
KAPPABRIDGES

AGICO

ADVANCED GEOSCIENCE INSTRUMENTS COMPANY

Supplement for KLY5-A and MFK1-FA Kappabridges for measurements of high-temperature variations of magnetic susceptibility.

General Description

The **CS4 High Temperature Furnace Apparatus** is optional attachment for **KLY5-A** and **MFK1-FA/A** Kappabridges designed for measurement of the temperature variations of low-field magnetic susceptibility of minerals, rocks and synthetic materials in the temperature range from **ambient temperature to 700°C**. Measurements can be performed under the protective argon atmosphere to prevent oxidation of measured specimen.

The apparatus consists of a non-magnetic furnace with a special platinum temperature sensor, a temperature control unit, and a cooling water reservoir. The specimen (up to 0.25 cm³ in volume) is placed in a silica glass vessel, heated by a platinum wire, and the temperature is measured by the temperature sensor.

The quasicontinuous measurement process is fully automated, being controlled by the software **Su-fyfe**. Data processing software **Cureval** serves for advanced analysis of thermomagnetic curves such as empty furnace measurement subtraction, Curie temperature estimation and separation of ferromagnetic and paramagnetic part of susceptibility.

CS4 Furnace Apparatus Comprising

- CS4 Temperature Control Unit
- Furnace
- Water Cooling Reservoir
- Temperature Sensor
- Argon Gas Flowmeter
- Power Supply Unit
- Specimen Vessels
- Set of Interconnecting Cables
- Measuring Software, CUREVAL Software
- User's Manual

Main Features

- Lowest detectable susceptibility change 1×10^{-7} SI.
- Measurement at high temperatures up to 700°C.
- Software controlled heating and cooling modes.
- Measurement in the air or argon atmosphere.
- 500 to 700 pairs of susceptibility and temperature measurements to define a thermomagnetic curve.

Technical specifications

Nominal specimen volume:	0.25 cm ³
Inner diameter of measuring vessel:	6 mm
Temperature range:	room temp up to 700°C
Accuracy of temperature determination:	$\pm 2^\circ\text{C}$
Sensitivity to susceptibility changes:	1×10^{-7} SI
Power requirements:	100 - 240 V, 50/60 Hz, 700 VA



LABORATORY INSTRUMENTS FOR MEASUREMENT OF MAGNETIC PROPERTIES OF ROCKS